

This form contains administration information for the MRI Control Room and Computer.

1 Keybox Code

6875 (last digits of system ID)

2 Kiosk Mode

Toggle in/out of kiosk mode via: Tab + Del + ‘Send To’

Credentials for XA software:

1. U: medadmin P: adm\$pwd\$4\$med\$
2. U: clinadmin P: clin\$pwd\$4\$med\$
3. U: medclient P: client\$pwd\$4\$med\$

3 Windows Mode

Ctrl + Esc

4 Reset Password

1. Log in as medadmin
2. Open ‘myExam Cockpit’
3. Select ‘Password Protection’

Old passwords will be loaded automatically, and it will be possible to set a new password.

5 Add DVS

Add DVS to appropriate location:

1. Enter Windows mode
2. K:
3. CustomerSeq
4. DiffusionVectorSets

Reference DVS in protocol:

1. Protocol Card
2. Diff
3. Diffusion Mode
4. Free